

## G1.01.Requirement\_profile

**Document:** G1.01.Requirement\_profile

**Release Date:** 05.07.2022

**Owner:** Head of R&D

**Template Version:** 6

# Table of contents

1. Summary	1
1.1. Responsible	1
1.2. Identifier	1
1.3. Description	1
1.3.1. What should be improved/solved with the requirement?	1
1.3.2. Value	1
1.3.3. Supported Devices	2
2. Use Cases	3
2.1. Use Case ARC-Req2_2-Uc1 (AF)	3
2.1.1. Risks	3
2.1.1.1. Risk ARC-Req2_2-Uc1-Risk1	3
2.1.1.1.1. Hazard	3
2.1.1.1.2. Harm	3
2.1.1.2. Risk ARC-Req2_2-Uc1-Risk2	3
2.1.1.2.1. Hazard	3
2.1.1.2.2. Harm	3
2.1.1.3. Risk ARC-Req2_2-Uc1-Risk3	3
2.1.1.3.1. Hazard	3
2.1.1.3.2. Harm	3
2.1.1.4. Risk ARC-Req2_2-Uc1-Risk4	4
2.1.1.4.1. Hazard	4
2.1.1.4.2. Harm	4
2.1.1.5. Risk ARC-Req2_2-Uc1-Risk5	4
2.1.1.5.1. Hazard	4
2.1.1.5.2. Harm	4
2.2. Use Case ARC-Req2_2-Uc2 (AFL)	4
2.2.1. Risks	4
2.2.1.1. Risk ARC-Req2_2-Uc2-Risk1	4
2.2.1.1.1. Hazard	4
2.2.1.1.2. Harm	4
2.2.1.2. Risk ARC-Req2_2-Uc2-Risk2	5
2.2.1.2.1. Hazard	5
2.2.1.2.2. Harm	5
2.2.1.3. Risk ARC-Req2_2-Uc2-Risk3	5
2.2.1.3.1. Hazard	5
2.2.1.3.2. Harm	5
2.2.1.4. Risk ARC-Req2_2-Uc2-Risk4	5
2.2.1.4.1. Hazard	5
2.2.1.4.2. Harm	5

2.2.1.5. Risk ARC-Req2_2-Uc2-Risk5 . . . . .	5
2.2.1.5.1. Hazard . . . . .	5
2.2.1.5.2. Harm . . . . .	5
2.3. Use Case ARC-Req2_2-Uc3 (sinus) . . . . .	5
2.3.1. Risks . . . . .	6
2.3.1.1. Risk ARC-Req2_2-Uc3-Risk1 . . . . .	6
2.3.1.1.1. Hazard . . . . .	6
2.3.1.1.2. Harm . . . . .	6
2.3.1.2. Risk ARC-Req2_2-Uc3-Risk2 . . . . .	6
2.3.1.2.1. Hazard . . . . .	6
2.3.1.2.2. Harm . . . . .	6
2.3.1.3. Risk ARC-Req2_2-Uc3-Risk3 . . . . .	6
2.3.1.3.1. Hazard . . . . .	6
2.3.1.3.2. Harm . . . . .	6
2.3.1.4. Risk ARC-Req2_2-Uc3-Risk4 . . . . .	6
2.3.1.4.1. Hazard . . . . .	6
2.3.1.4.2. Harm . . . . .	6
2.4. Use Case ARC-Req2_2-Uc4 (noise) . . . . .	7
2.4.1. Risks . . . . .	7
2.4.1.1. Risk ARC-Req2_2-Uc4-Risk1 . . . . .	7
2.4.1.1.1. Hazard . . . . .	7
2.4.1.1.2. Harm . . . . .	7
2.5. Use Case ARC-Req2_2-Uc5 (no supported ecg format) . . . . .	7
2.5.1. Risks . . . . .	7
2.5.1.1. Risk ARC-Req2_2-Uc5-Risk1 . . . . .	8
2.5.1.1.1. Hazard . . . . .	8
2.5.1.1.2. Harm . . . . .	8

# 1. Summary

## 1.1. Responsible

**Jens Meder** is the contact person for this requirement. The contact person accompanies this requirement over the entire life cycle.

## 1.2. Identifier

ARC-Req2\_2 : Analyse ecg data

## 1.3. Description

### 1.3.1. What should be improved/solved with the requirement?

- Analyze an ecg including all available ecg channels
  - Support for multiple [ECG devices](#)
  - Detect and annotate episodes of Atrial Fibrillation (AF) and Atrial Flutter (AFI) of 30s and more
  - Detect and annotate signal noise and artefacts in each channel
  - Calculate the atrial fibrillation burden (AF Burden) in terms of
    - the longest episode of AF/AFI measured in seconds
    - the percentage of AF/AFI, e.g., 14% for a total of 42s of af in a 300s recording
  - Calculate the Heart Rate Variability (HRV) using root mean square of successive differences (RMSSD) in milliseconds
  - Calculate the signal quality in terms of net time without artefacts or signal noise
  - Calculate the heart rate (exclude noisy areas to avoid measurement errors)
  - Analysis result is uniquely identifiable (without referencing the specific patient)

### 1.3.2. Value

- Analysis result include following information:
  - AF was correctly detected and annotated
  - Noise was correctly detected and annotated
  - Calculated AF burden
  - Calculated HRV
  - Calculated signal quality
  - Calculated heart rate
- Authorized user can access the correct result
- Authorized user generally do not need more than 20 seconds to access the result

- The user needs maximum 3 seconds to classify the result.
- A printable document format is available.
- Notes for proposal:
  - The Analysis will have a medical purpose. This makes the following steps necessary:
    - More details will be part of the validation plan and clinical evaluation
    - Precise definition is therefore made in the corresponding process activities

### 1.3.3. Supported Devices

- Bittium Faros

## 2. Use Cases

### 2.1. Use Case ARC-Req2\_2-Uc1 (AF)

- Personas
  - Employee in a private medical practice
  - Desktop PC of the practice
- Preconditions
  - ECG data of a patient with AF (atrial fibrillation) is recorded and available for analysis
- Procedure
  - The ECG data is provided for analysis.
- Expected result
  - The employee interprets the visualized analysis result as AF.

#### 2.1.1. Risks

##### 2.1.1.1. Risk ARC-Req2\_2-Uc1-Risk1

###### 2.1.1.1.1. Hazard

- False positive analysis result

###### 2.1.1.1.2. Harm

- Further treatment due to false positive analysis is performed

##### 2.1.1.2. Risk ARC-Req2\_2-Uc1-Risk2

###### 2.1.1.2.1. Hazard

- Incorrect assignment of the analysis report

###### 2.1.1.2.2. Harm

- A patient with AF is not recognized

##### 2.1.1.3. Risk ARC-Req2\_2-Uc1-Risk3

###### 2.1.1.3.1. Hazard

- Presentation of a non-representative 5-minute episode

###### 2.1.1.3.2. Harm

- A patient with AF is not recognized and will not be treated accordingly

#### 2.1.1.4. Risk ARC-Req2\_2-Uc1-Risk4

##### 2.1.1.4.1. Hazard

- Wrong interpretation of the analysis result due to untrained personnel

##### 2.1.1.4.2. Harm

- A patient with AF is not recognized and will not be treated accordingly

#### 2.1.1.5. Risk ARC-Req2\_2-Uc1-Risk5

##### 2.1.1.5.1. Hazard

- Wrong interpretation of the analysis result due to misleading presentation or unreadable analysis result

##### 2.1.1.5.2. Harm

- A patient with AF is not recognized and will not be treated accordingly

### 2.2. Use Case ARC-Req2\_2-Uc2 (AFL)

- Personas
  - Employee in a private medical practice
  - Desktop PC of the practice
- Preconditions
  - ECG data of a patient with AFL (atrial flutter) is recorded and available for analysis
- Procedure
  - The ECG data is provided for analysis.
- Expected result
  - The employee interprets the visualized analysis result as AFL.

#### 2.2.1. Risks

##### 2.2.1.1. Risk ARC-Req2\_2-Uc2-Risk1

###### 2.2.1.1.1. Hazard

- False positive analysis result

###### 2.2.1.1.2. Harm

- Further treatment due to false positive analysis is performed

#### **2.2.1.2. Risk ARC-Req2\_2-Uc2-Risk2**

##### **2.2.1.2.1. Hazard**

- Incorrect assignment of the analysis report

##### **2.2.1.2.2. Harm**

- A patient with AFL is not recognized

#### **2.2.1.3. Risk ARC-Req2\_2-Uc2-Risk3**

##### **2.2.1.3.1. Hazard**

- Presentation of a non-representative 5-minute episode

##### **2.2.1.3.2. Harm**

- A patient with AFL is not recognized and will not be treated accordingly

#### **2.2.1.4. Risk ARC-Req2\_2-Uc2-Risk4**

##### **2.2.1.4.1. Hazard**

- Wrong interpretation of the analysis result due to untrained personnel

##### **2.2.1.4.2. Harm**

- A patient with AFL is not recognized and will not be treated accordingly

#### **2.2.1.5. Risk ARC-Req2\_2-Uc2-Risk5**

##### **2.2.1.5.1. Hazard**

- Wrong interpretation of the analysis result due to misleading presentation or unreadable analysis result

##### **2.2.1.5.2. Harm**

- A patient with AFL is not recognized and will not be treated accordingly

### **2.3. Use Case ARC-Req2\_2-Uc3 (sinus)**

- Personas
  - Employee in a private medical practice
  - Desktop PC of the practice
- Preconditions
  - ECG data of a patient with a normal sinus rhythm is recorded and available for analysis
- Procedure



- The ECG data is provided for analysis.
- Expected result
  - The employee interprets the visualized analysis result as a normal ecg without arrhythmia.

## **2.3.1. Risks**

### **2.3.1.1. Risk ARC-Req2\_2-Uc3-Risk1**

#### **2.3.1.1.1. Hazard**

- False negative analysis result

#### **2.3.1.1.2. Harm**

- Further treatment due to false negative analysis is performed

### **2.3.1.2. Risk ARC-Req2\_2-Uc3-Risk2**

#### **2.3.1.2.1. Hazard**

- Incorrect assignment of the analysis report

#### **2.3.1.2.2. Harm**

- A patient with sinus ecg rhythm is not recognized as such, and will be treated inappropriately.

### **2.3.1.3. Risk ARC-Req2\_2-Uc3-Risk3**

#### **2.3.1.3.1. Hazard**

- Wrong interpretation of the analysis result due to untrained personnel

#### **2.3.1.3.2. Harm**

- A patient with sinus ecg rhythm is not recognized as such, and will be treated inappropriately.

### **2.3.1.4. Risk ARC-Req2\_2-Uc3-Risk4**

#### **2.3.1.4.1. Hazard**

- Wrong interpretation of the analysis result due to misleading presentation or unreadable analysis result

#### **2.3.1.4.2. Harm**

- A patient with sinus ecg rhythm is not recognized as such, and will be treated inappropriately.

## 2.4. Use Case ARC-Req2\_2-Uc4 (noise)

- Personas
  - Employee in a private medical practice
  - Desktop PC of the practice
- Preconditions
  - ECG data with noise of a patient is recorded and available for analysis
- Procedure
  - The ECG data is provided for analysis.
- Expected result
  - The employee interprets the visualized analysis result as a ecg with noise.

### 2.4.1. Risks

#### 2.4.1.1. Risk ARC-Req2\_2-Uc4-Risk1

##### 2.4.1.1.1. Hazard

- Incorrect assignment of the analysis report

##### 2.4.1.1.2. Harm

- A patient with AF/AFL is not recognized

## 2.5. Use Case ARC-Req2\_2-Uc5 (no supported ecg format)

- Personas
  - Employee in a private medical practice
  - Desktop PC of the practice
- Preconditions
  - ECG data of a patient is recorded and available for analysis, however the ecg data is not supported.
- Procedure
  - The ECG data is provided for analysis.
- Expected result
  - The employee interprets the visualized analysis result as a ecg that is not supported and cannot be analysed.

### 2.5.1. Risks

2.5.1.1. Risk ARC-Req2\_2-Uc5-Risk1

2.5.1.1.1. Hazard

No assignment of ecg data

2.5.1.1.2. Harm

A patient with AF/AFL is not recognized

Approval

Approver	Date and signature
PO	
R&D	